



..... Testing • Research • Consulting

Preset Values:

- ➔ Materials: Vacuum-resistant, metallic and non-metallic solids and powders. Magnetic materials can also be analyzed.
- ➔ Dimension: Sample diameter of max. 100 mm; height of max. 20 mm, larger samples can be reduced to the maximum dimensions.
- ➔ Sample Handling: Ideally, the samples are not touched by hand and packed in common aluminium foil for the transportation.

Contacts

Dr. Roman Heuberger
 Manager Surfaces & Tribology
 Phone +41(0)32 644 14 74
roman.heuberger@rms-foundation.ch

Olivier Loeffel
 XPS & SEM Operator
 Phone +41(0)32 644 14 55
olivier.loeffel@rms-foundation.ch

Please discuss with us! We will be happy to advise you.
Or ask for our service catalogue. You will find this and other information on our website as well.



S
T
S
SCHWEIZERISCHER PRUFSTELLENDIENST
 SERVICE SUISSE D'ESSAIS
 SERVIZIO DI PROVA IN SVIZZERA
 SWISS TESTING SERVICE

ISO 9001
 BUREAU VERITAS
 Certification
 N° CH1099316



..... Testing • Research • Consulting

XPS @ RMS



X-Ray Photoelectron Spectroscopy (XPS)

Chemical surface analysis for

- ➔ Cleanliness analysis
- ➔ Analysis of contaminants and residues
- ➔ Monitoring of surface modifications

Dr. h. c. Robert Mathys Stiftung • RMS Foundation
 Bismattstrasse 12 • P.O. Box 203 • CH-2544 Bettlach
 Phone +41(0)32 644 14 00 • Fax +41(0)32 644 11 76 • www.rms-foundation.ch

How does XPS work?

In an X-ray photoelectron spectrometer electrons are excited with X-rays so strongly that they leave their atom and eventually the sample surface as well. The energy of these photoelectrons is analysed and hence their binding energy is calculated. This permits to determine quantitatively the chemical composition in the top 5-10 nm of the surface.

This technique detects all elements except hydrogen and helium and their binding states. The detection limit is approximately 0.1 at%, which corresponds to about 1 ng/cm² on a surface.




Figure: Detection principle.

Applications

XPS is a non-destructive technique to analyse the

- Cleanliness of surfaces
- Chemical composition of residues and spots
- Chemical modification after surface treatments (e.g. for research)

Our costumers are mostly from the medical technology field but also from other fields such as watch, coating, electro- & semi conductor industry.

Equipment

Axis Nova from Kratos Analytical, Manchester, UK.

Our spectrometer is equipped with a monochromatic AlK α X-ray source, charge neutraliser, argon ion sputtering, a hemispherical analyser and a 2D-detector that allows for the acquisition of chemical maps with a lateral resolution of 3 μ m, which is the strength of our instrument!

Types of Analysis & Costs

XPS Cleanliness Analysis / Standard Surface Analysis

Such an analysis includes the measurement of survey spectra, data processing and a report (German or English). The report consists of the spectra, a table with the quantitative chemical composition and information about the chemical bindings (e.g. distinction between sulphates and sulphides).

Table: Costs for XPS Cleanliness Analysis / Standard Surface Analysis

No. of analysed spots	Costs [CHF] *
1	550
2	770
3	930
4	1070
5	1210
6	1350
7	1510
8	1680
9	1840
10	2000
further Spots	+200/Spot

* XPS Analyses of special spots, complicated or magnetic samples need more time to analyse and thus the additional expenses are charged (approximately 10-50% higher than for standard analyses).

Special XPS Investigations

- Detailed spectra: Acquisition of detailed spectra of selected elements for a more accurate chemical analysis.
- Angle resolved XPS: Analysis of tilted samples in order to determine the depth distribution of the elements in the top 10 nm.
- Depth profiling: Applying argon ion sputtering, the depth profile between 10 nm and 1 μ m depth can be determined.
- Imaging XPS: Chemical map with the distribution of elements or oxidation states with a resolution of 3 μ m!

For these special investigations the actual expenses are charged (use of XPS: 430 CHF/h, personnel cost for the measurement, analysis and report: 170-200 CHF/h). Please ask us for a quotation.

Accreditation & Norms

Our XPS services are accredited according to ISO 17025 and certified according to ISO 9001. The Norms ISO 10993-18 and ASTM F2847 recommend XPS for the in-situ analysis of the chemical surface composition of medical devices (for organics, inorganics and insoluble particulates).

- ISO 17025: General requirements for the competence of testing and calibration laboratories
- ISO 9001: Quality management systems - Requirements
- ISO 10993-18: Biological evaluation of medical devices - Part 18: Chemical characterization of materials
- ASTM F2847: Standard Practice for Reporting and Assessment of Residues on Single Use Implants