

Protection of electronics for reprocessible surgical devices

A. Lamibrac, G. Schaffner

Pre-development unit, Turck duotec SA, CH-2800 Delémont, Switzerland;
adrien.lamibrac@turck.com, gregory.schaffner@turck.com



A complete stand alone solution with connections like PTFE wires or gold connector up to over 1000 cleaning and steam sterilization cycles for your surgical devices

Background & Goals

- ❖ Miniaturized electronics like sensors is increasingly used in surgical devices
- ❖ Reprocessing (cleaning + steam sterilization) of surgical device is harsh for electronics and can cause early failure
- ❖ Building a housing around the electronics and providing an efficient barrier effect against reprocessing is necessary if the device has to be reprocessed several times
- ❖ An electrical interface (i.e. PTFE wires or gold pins) through the housing is a potential “weakest link” in the barrier effect
- ❖ Understanding and managing the failure mechanism of these interface is critical to improve the resistance of electronics against reprocessing
 - ➔ Use of accelerated ageing to simulate reprocessing cycles
 - ➔ IP testing to detect any failure of the interface

Experimental



Figure 1: View of the Device Under test (DUT) used in this study

- ❖ Device Under Test (DUT)'s construction
 - A test PCB is equipped with a golden pin perpendicularly soldered on top of it
 - Through adequate pre-treatment, gold's surface is prepared for enhanced bonding properties when overmolded
 - The assembly is overmolded with our specific compound for reprocessible devices
- ❖ The DUT is going through accelerated ageing in a two chamber oven from 0 to 150° C. The extended temperature range compensate the lower heat transfer rate compared to sterilization.
- ❖ Interfaces between different material used in the device construction will be periodically loaded and induce mechanical fatigue of the bond between them.
- ❖ The DUT is then immersed under 1 meter of water for 30 minutes
- ❖ A dielectric test in water under 1500 V is performed to measure if any water ingress happened during the immersion



Figure 2: View of the two chamber oven where the DUT are going through accelerated ageing 0 - 150°C




Figure 3: Scheme of the dielectric test and potential ingress point inducing current leakage

Results

- ❖ Theoretical work results
 - Potting or overmolding **removes microclimate** but **adds mechanical stress** on the PCB and the electronic components
 - To avoid mechanical failure of the solder joints or PCB, the **balance between mechanical properties of the polymer and max shear strain of the soldering joints** of each components has to be found
 - During **steam sterilization**, steam condensates onto the device's surface. Water has a **heat transfer rate 30 times higher** than steam.
 - Temperature gradient throughout the sterilized device are important causing **efforts and mechanical fatigue**, especially at the **interfaces between two material**
- ❖ Practical results
 - The initial current leakage of the DUT during dielectric test after an IP67 immersion is around 0.010 mA
 - After pre-ageing the DUT up to 300 shock cycles or sterilization, **dielectric test results stay low and steady around 0.010 mA**





Figure 4: Leakage current measurement after different ageing processes
No significant rise of leakage current indicates a stable interface quality

Conclusions

- ❖ Electronic can be **reliably protected** from reprocessing
- ❖ Depending on the electrical interfaces needed, the performance can go up to **over 1000 reprocessing cycles**
- ❖ Metallic pins can be used as connection elements while keeping the **necessary barrier effect** against reprocessing
- ❖ The Stand-alone protection for electronic concept allows more **design freedom for the designer of the complete surgical device**