

Specialized Calcium Phosphate Materials for Research Applications

RMS Foundation has been active in the research of calcium phosphate bioceramics for more than 35 years. During this time, our experts have developed materials with outstanding physical and chemical properties. As a testing, research, and consulting institute, we offer our expertise to partners from academia and industry. Our services include transfer and implementation of specialized manufacturing processes to the clients' facilities, consulting to improve existing processes, and supply of raw materials and samples for research projects.

High-purity α -TCP powder

Our in-house produced α -TCP (α -Ca $_3$ (PO $_4$) $_2$) reaches a phase purity of 98% or more. The material can be mixed with water to form a cement paste. A passivated variation is more convenient to use in production processes requiring longer processing times, e.g. for 3D printing / robocasting. A standard production batch yields up to 3 kg, and the process allows fine adjustments to the Ca:P ratio and the final particle size.

Amorphous CaP / nano-α-TCP

Amorphous calcium phosphate (ACP) and nanoparticulate α -TCP powders of high chemical purity are produced by precipitation synthesis and thermal treatment. When used in pastes for cement applications or 3D printing / robocasting, these nanoparticulate materials allow to optimize the flow, injectability, and setting characteristics. ACP reacts quickly with aqueous solutions, whereas the reaction of nano- α -TCP is initially delayed, but then proceeds more rapidly than our regular high-purity α -TCP.

Ultra-pure β-TCP

Ultra-pure β -TCP (β -Ca $_3$ (PO $_4$) $_2$) is produced from highly pure raw materials in an optimized precipitation process that allows us to control the final phase purity and to incorporate specific chemical substitutions [1]. Due to the precise control over the chemical composition, this material can be used to investigate the effects of ionic substitutions


on biological, chemical, and physical properties of the material without interference from impurities.


Dense β-TCP

Using our ultra-pure β -TCP and an optimized slip-casting process, we can produce disks of highly densified β -TCP. Other phase compositions (HA, BCP) can be produced upon request. Due to a density of >99%, polished surfaces are extremely smooth and allow to visualize osteoclast pits by electron or confocal microscopy. The size of our disks is matched to 96 well plates, which makes them an ideal substrate for cell experiments. Modifications to provide calcium or phosphate enriched surface compositions are also possible [2].

Material characterizations

RMS Foundation offers a comprehensive suite of physical and chemical analytical services for characterizations of calcium phosphate materials and products according to the relevant ISO and ASTM standards.





Highly densified grain structure of a β -TCP disk before (top) and after (bottom) partial resorption by macrophages [3].

[1] Le Gars Santoni, B. et al., J. Eur. Ceram. Soc. 41 (2021), 1683-1694.

[2] Döbelin, N. et al., J. Eur. Ceram. Soc. 40 (2020), 6095-6106.

[3] Sblendorio, G. A., PhD thesis, EPFL, Switzerland 120211.

Newsletter 38

α-TCP:

- Phase purity: >98 wt-%
- Impurities of β-TCP, β-CPP, HA: 0
 1 wt-%
- Batch size: 2.5 3 kg

Nano-α-TCP:

- Phase purity: >98 wt-%
- Crystallite size: ~80 nm
- Specific surface area: 12 15 m²/g
- Chemical impurities: <20 ppm
- Batch size: 15 20 g

Ultra-pure β-TCP:

- Phase purity: >99%
- Chemical impurities: <20 ppm
- Substitutions: Sr, Mg, Cu, Gd...
- Batch size: 15 20 g

Dense β-TCP disks:

- Phase purity: >99%
- Chemical impurities: <20 ppm
- Density: >99%
- Ø ~5 mm, height ~2 mm

Discuss your questions with us! We will be happy to advise you.

Contact for more information:

Dr. Nicola Döbelin Phone +41 32 644 20 41 nicola.doebelin@rms-foundation.ch

More information as well as our service catalogue can be found on our website.

The RMS Foundation has been certified according to ISO 9001 and is an ISO/IEC 17025 (type C) accredited laboratory.

Subscribe to our mailing list and read more RMS newsletters covering other topics.

